

## OsseoSpeed™ TX Profile

After a tooth extraction or tooth loss, alveolar ridge remodeling occurs as the bone heals. The remodeling has been shown to be more pronounced on the buccal side than on the lingual, often resulting in a sloped ridge in a lingual to buccal direction<sup>1,2</sup>. When placing a regular implant in a healed sloped ridge the implant can either be installed in level with the lingual marginal bone or in level with the buccal marginal bone. However, this will result in exposed threads on the buccal aspect or unsupported lingual bone, both situations being non-optimal from a long-term esthetical point of view. A more favourable solution in a sloped ridge situation would be to use an implant with a sloped shoulder that fits the bony anatomy, optimizing the bone support thus reducing the risk of losing unsupported marginal bone. The OsseoSpeed™ TX Profile, an implant with a sloped shoulder, was developed in order to facilitate implant treatment in sloped ridge situations.

A clinical multicenter study evaluating OsseoSpeed™ Profile implants in single-tooth positions in patients with healed sloped ridges is ongoing. Sixteen-week re-entry data, measuring buccal and lingual bone levels, showed well maintained marginal bone levels<sup>3</sup>. The patients have been followed for 1 year, no implants have been lost, and the results indicate that the OsseoSpeed Profile implant is a predictable treatment option in cases where the alveolar crest is sloped in a lingual to buccal direction<sup>4,5</sup>.

An implant can also be placed into the fresh extraction socket. However, studies have shown that despite an implant in place, similar remodeling of the alveolar crest occurs, resulting in a sloped ridge in a lingual to buccal direction<sup>6</sup>. Placement of OsseoSpeed™ Profile implants using an immediate implant placement protocol, have shown a nice correlation to the sloped ridge after remodeling and thereafter maintained circumferential marginal bone levels and good esthetic results after 1 year<sup>7</sup>.

## References

Reprints can be ordered from references marked with ID No.  
To read more Scientific Reviews please see: [www.astratechdental.com](http://www.astratechdental.com)

1. Araujo MG, Lindhe J. Dimensional ridge alterations following tooth extraction. An experimental study in the dog. J Clin Periodontol 2005;32(2):212-8. [Abstract in PubMed](#)

2. Van der Weijden F, Dell'Acqua F, Slot DE. Alveolar bone dimensional changes of post-extraction sockets in humans: a systematic review. J Clin Periodontol 2009;36(12):1048-58. [Abstract in PubMed](#)

3. Noelken R, Donati M, Fiorellini J, Gellrich. Maintenance of marginal bone support at OsseoSpeed™ Profile implants. Clin Oral Impl Res 2010;21(10):#201:1079. [Abstract in PubMed](#)

4. Fiorellini J, Parker W, Nölken R, et al. Marginal bone support at OsseoSpeed™ Profile implants. Academy of Osseointegration, Washington, DC, 2011:#126 2011.

5. Noelken R. Marginal bone maintenance and soft tissue esthetics with the Astra Tech OsseoSpeed™ Profile implant: a case report. Dentale Implantologie & Parodontologie 2011;March.

6. Sanz M, Cecchinato D, Ferrus J, et al. A prospective, randomized-controlled clinical trial to evaluate bone preservation using implants with different geometry placed into extraction sockets in the maxilla. Clin Oral Implants Res 2010;21(1):13-21. [Abstract in PubMed](#)

7. Noelken R, Neffe B, Kunkel M, Wagner W. Maintenance of marginal hard and soft tissue support at immediately provisionalized OsseoSpeed Profile implants – 1-year results. Clin Oral Implants Res 2011;22(9):#180:959. [Abstract in PubMed](#)

